



## Testimonials

*"Sentia is so easy. I save a bunch of time, a bunch of money, I do less work and it's perfectly accurate. Every winemaker needs one, there's no question."*

**Sandor Faludi, Winemaker, Prohibition Cellars  
Woodinville, WA, USA**

*"For all the tests this unit will do makes it a smart, affordable choice for any size winery..."*

**Chris Eberle, Winemaker, Eberle Winery, Paso Robles, CA, USA**

*"I hardly need lab equipment anymore with my new Sentia. Faster, easier, and more accurate. This thing is the future."*

**Tyler Farnsworth, Winemaker, Laterus Winery  
Snohomish, WA, USA**

*"For us, this innovative and simple to use device has taken our task of free SO<sub>2</sub> monitoring from a two-week process down to two days."*

**Don Schroeder, Winemaker, Sea Smoke  
Santa Rita Hills, Santa Barbara, CA, USA**

*"The Sentia device has been a game-changer for me... I can check in on individual barrels or tanks or do composite checks any time that I'd like at a fraction of the cost of having a third party lab run the same analysis."*

**Larry Schaffer, Owner & Winemaker, Tercero Wines  
Los Olivos, CA, USA**

*"I cannot even begin to calculate the amount of hours saved in the lab because of the Sentia."*

**Riley Hollenbach, Winemaker, Bartier Bros. Vineyard & Winery  
Oliver, British Columbia, Canada**

*"Using the Sentia has been fantastic! It has meant that I am able to assess our grape juice and wine more frequently as well as shortened the turn around time with results, ultimately saving on time and enabling better wine making decisions."*

**Brad Atkins, Owner/Vigneron, Steels Gate Wines  
Dixons Creek, VIC, Australia**

*"My father-in-law has been making wine for 40 years and has always been sceptical of new technology. When he used Sentia the first few times he was instantly converted!"*

**Stephanie Helm & Ben Osborne, Winemaker & Viticulturalist,  
Vintners Daughter - Murrumbateman, NSW, Australia**



<https://sentiaanalysis.com/>

**For a free demonstration call:**

Designed by

Universal Biosensors Pty Ltd

1 Corporate Avenue, Rowville, 3178, Victoria, Australia



[www.universalbiosensors.com](http://www.universalbiosensors.com)

[SentiaSales@universalbiosensors.com](mailto:SentiaSales@universalbiosensors.com)



**Measure 6 key wine analytes in minutes**

Now available in North America, Australia, New Zealand, Europe, Chile and South Africa



**SENTIA™**

**The most convenient wine analysis tool available**



**Portable**  
test in the winery



**Touch screen**  
for ease of use



**Easy, intuitive**  
no special training  
required



**Small sample size**  
less wastage



**Quick results**  
fast decisions



**Auto-calibrating**  
peace of mind



**Affordable**  
save on costs

## Product range

<b>Meter</b>	Product code: 90200
<b>Interface</b>	Wireless (2.4 GHz)
<b>Automatic power off</b>	Configurable: disabled; 10 mins, 20 mins, 30 mins
<b>Dimensions</b>	147 x 84 x 32 mm (5.79 x 3.31 x 1.26 in)
<b>Weight</b>	210 g (7.41 oz)
<b>Meter Case (+stylus)</b>	Product code: 90210
<b>Dimensions</b>	161.5 x 92 x 41.5 mm (6.36 x 3.62 x 1.63 in)
<b>Weight</b>	89 g (3.14 oz)
<b>Acetic Acid</b>	Product code: 30730
<b>Sample type</b>	Red or white wine
<b>Measurement range</b>	0.1 to 1.5 g/L
<b>Test time</b>	<3.5 min
<b>Dilution required</b>	Yes (1 part sample, 3 parts deionized/ distilled water)
<b>Free SO<sub>2</sub></b>	Product code: 30230
<b>Sample type</b>	Post-fermentation red or white wine
<b>Measurement range</b>	3 to 50 mg/L
<b>Test time</b>	<1 min
<b>Dilution required</b>	No
<b>Fructose</b>	Product codes: 30530 (strips); 91004 (diluent)
<b>Sample type</b>	Red or white wine
<b>Measurement range</b>	0.1 to 10 g/L
<b>Test time</b>	<2 min
<b>Dilution required</b>	Yes (1 part sample, 4 parts Sentia fructose buffer)
<b>Glucose</b>	Product code: 30330
<b>Sample type</b>	Red or white wine
<b>Measurement range</b>	0.1 to 10 g/L
<b>Test time</b>	<1 min
<b>Dilution required</b>	No
<b>Malic Acid</b>	Product codes: 30430 (strips); 91002 (diluent)
<b>Sample type</b>	Red or white wine
<b>Measurement range</b>	0.05 to 5 g/L
<b>Test time</b>	<1 min
<b>Dilution required</b>	Yes (1 part sample, 4 parts Sentia malic acid buffer)
<b>Titrateable Acidity (pH 7.00 &amp; 8.20)</b>	Product code: 30630
<b>Sample type</b>	Red or white wine
<b>Measurement range</b>	3 to 10 g/L
<b>Test time</b>	<1 min
<b>Dilution required</b>	No



Easy-to-follow touch-screen prompts



Test in the winery



Small sample required > 8 µL to 100 µL



Results detected within minutes